

Distributed By

Medical Refrigeration .Co.Uk

JustCatering.com Ltd
 10 Jessop Avenue, Norwood Green,
 Middlesex, UB2 5UY, UK
 Tel & Fax: +44 (0)8450034271
 Sales Direct Line: 07006342365
 E-mail: sales@medicalrefrigeration.co.uk
 Web: http://www.medicalrefrigeration.co.uk

First Aid Cabinets

Sun-FAC1	Small First Aid Cabinet
Sun-FAC3LW	Medium Size First Aid Cabinet with Lock in - White
Sun-FAC3LS	Medium First Aid Cabinet with Lock in - Silver
Sun-FAC4W	Large First Aid Cabinet with Lock in - White
Sun-FAC4S	Large First Aid Cabinet with Lock in - Silver

39



SUNFLOWER Medical Furniture Contemporary by Design

Bins

New




Our range of hands free, flame retardant and clinical waste bins are designed for the multi-function disposal of waste. Usage is designated by an appropriate coloured label (blue for linen/laundry and black for domestic) and colour coded yellow lid and pedal label for clinical waste versions.

Sizes



Waste Identification



Product Features

hands free

Hygienic Operation

coloured pedal labels

Identifies Bin Usage/Contents



Silent, Dampened Closing Lid

3 sizes

tested to HTM83 standard



- Three sizes, 20, 50 and 70 litre capacity
- Hands free operation, foot pedal must be used.
- Corrosion resistant galvanised bodies with stainless steel foot pedal and mechanism
- Sealed base prevents leakages and facilitates easy cleaning
- Integrated wall mounting bracket
- Elasticated cord hold bags firmly in place, while the nylon loop makes changing bags easier
- None clinical bins supplied with choice of foot pedal labels
- Clinical bins have yellow lids and clinical waste foot pedal label
- Designed and tested to Healthcare Technical Memorandum (HTM83)

Dimensions

20Ltr	50Ltr	70Ltr
51cm(H) x 32cm(W) x 33cm(D)	59cm(H) x 42cm(W) x 43cm(D)	78cm(H) x 42cm(W) x 43cm(D)

40