


product information

 Electric or Hydraulic height adjustment 45-96cm	 Maximum user weight (35st.)	 Gas strut assisted multi-position head and foot sections	 Back rest and leg rest adjustment +70° / -35°	 Easy to use handset (Electric models only)
 Retractable positioning castors	 Adjustable height levelling foot	 Fully upholstered head & foot sections	 Double layer density foam for support & comfort	 Luxury 65mm deep contoured upholstery
 Anti-bacterial/anti-microbial vinyl	 Paper roll holder (supplied as standard)	 Easy clean design and materials	 3 year warranty on all moving components	 Lifetime warranty on frame

accessories \* these accessories can only be fitted when ordered with the couch















 <b>Sun-FVHC/LEG1</b> Leg Supports	 <b>Sun-FVHC/PHB1</b> Phlebotomy Arm	 <b>Sun-FVHC/FSW1*</b> Foot Switch	 <b>Sun-FVHC/SUP1</b> Gynae Stirrups	 <b>Sun-CHP1</b> Head Pillow
 <b>Sun-FVHC/CBH1*</b> Breathing Hole	 <b>Sun-FVHC/CSR</b> Folding Side Rails (pair)	 <b>Sun-FVHC/IV2</b> IV Pole Attachment	 <b>Sun-FVHC/RAIL</b> Medi-Rail Attachment	 <b>Sun-FVHC/BB*</b> Battery Backup
 <b>Sun-FVHC/HEO*</b> Electrically Operated Head Section Upgrade	 <b>Sun-FVHC/FEO*</b> Electrically Operated Foot Section Upgrade	 <b>Sun-INST1*</b> Installation Service & Packaging Recycling	 <b>Sun-DISPOSAL1*</b> Removal & Disposal of Old Couch	

dimensions



anti-bacterial vinyl upholstery colours - please specify when ordering

please see inside back cover for full specification of our vinyl range


														
black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Fusion Classic - 3 Section Variable Height Couches

The new generation of Fusion Classic has been designed to assist in making treatment procedures even easier and more efficient than ever. Fusion Classic is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head & foot sections with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 225kg



**Optional**  
Electric Head & Foot Sections (electric models only)



Normally delivered within 7-10 working days (UK Only)



Code: **Sun-FVHC3E**  
Fusion Classic 3 Section Electric Couch (Shown with optional accessories)



Code: **Sun-FVHC3H**  
Fusion Classic 3 Section Hydraulic Couch (Shown with optional pillow)



Hydraulic height adjustment

product information

 Electric or Hydraulic height adjustment 45-96cm	 Maximum user weight (35st.)	 Gas strut assisted multi-position head section	 Back rest adjustment +70° / -35°	 Easy to use handset (Electric models only)
 Retractable positioning castors	 Adjustable height levelling foot	 Fully upholstered head section	 Double layer density foam for support & comfort	 Luxury 65mm deep contoured upholstery
 Anti-bacterial/anti-microbial vinyl	 Paper roll holder (supplied as standard)	 Easy clean design and materials	 3 year warranty on all moving components	 Lifetime warranty on frame

accessories \* these accessories can only be fitted when ordered with the couch

 <b>Sun-FVHC/LEG1</b> Leg Supports	 <b>Sun-FVHC/PHB1</b> Phlebotomy Arm	 <b>Sun-FVHC/FSW1*</b> Foot Switch	 <b>Sun-FVHC/SUP1</b> Gynae Stirrups	 <b>Sun-CHP1</b> Head Pillow
 <b>Sun-FVHC/CBH1*</b> Breathing Hole	 <b>Sun-FVHC/CSR</b> Folding Side Rails (pair)	 <b>Sun-FVHC/IV2</b> IV Pole Attachment	 <b>Sun-FVHC/RAIL</b> Medi-Rail Attachment	 <b>Sun-FVHC/BB*</b> Battery Backup
 <b>Sun-FVHC/HEO*</b> Electrically Operated Head Section Upgrade	 <b>Sun-INST1*</b> Installation Service & Packaging Recycling	 <b>Sun-DISPOSAL1*</b> Removal & Disposal of Old Couch		


dimensions



185cm(L) 65cm(W)



anti-bacterial vinyl upholstery colours - please specify when ordering  
please see inside back cover for full specification of our vinyl range


														
black	grey	beige	primrose	walnut	red wine	salmon	lavender	navy	mid blue	sky blue	cool blue	green	mint	white

Fusion Classic - 2 Section Variable Height Couches

The new generation of Fusion Classic has been designed to assist in making treatment procedures even easier and more efficient than ever. Fusion Classic is simply, strong, stable, aesthetically pleasing and easy to clean.



Fully upholstered head section with gas assisted adjustment



Electric or Hydraulic height adjustment



Increased maximum user weight to 225kg



**Optional**  
Electric Head Section  
(electric models only)



Normally delivered within 7-10 working days (UK Only)



Code: **Sun-FVHC2E**  
Fusion Classic 2 Section Electric Couch



Code: **Sun-FVHC2H**  
Fusion Classic 2 Section Hydraulic Couch



Hydraulic height adjustment