

Pharmacy Refrigerator

MED140



Stackable with kit
M4-UCRSTACK
£76.03



How to Set Alarm
Thermometer
Video



For convenience the thermometer sensor has been fitted to the cabinet ready for use. Supplied with 2x monthly Temperature log sheets ready for immediate use.




image: Freepik.com

In anticipation of approved vaccines for **Covid-19** and subsequent roll out of mass vaccination we at Pentland have developed the Blizzard Pharmacy refrigerator.

We believe the demand for pharmacy refrigeration will be huge.

The Blizzard Pharmacy refrigerator has been validated* independently to operate at **+2°C to +8°C** correctly at a very competitive cost.

The unit can be stacked (with stacking kit) to give double storage capacity on the same footprint.

Apart from Covid-19 related projects this unit is perfect for:

- Doctors surgeries
- Pharmacies
- Dentists
- Hospitals
- Veterinary surgeries

* Validated under specific operating conditions, site specific validation is recommended (NOT INCLUDED). Consult your validation technician / supplier.

FEATURES

- White laminated exterior; ABS interior
- Secondary alarm thermometer accuracy of +/- 0.5°C
- Audible and visual temperature alarm - HI/LO alarm
- MAX and MIN temperature monitoring
- Battery powered temperature monitoring in case of power failure
- Height adjustable from 819 – 850mm
- Stackable with stacking kit (M4-UCRSTACK)
- Reversible door
- Lock fitted as standard
- Adjustable shelves
- Replaceable door gasket

TECHNICAL SPECIFICATIONS

Controller:	Digital controller with temperature display. Pre-set and locked to operate at a set point of 5°C. Secondary alarm thermometer with fitted sealed glycol sensor.
Cooling:	Fan assisted
Defrost:	Automatic
Refrigerant:	R600a
Climate Class:	4
Temperature Range:	+2°C ~ +8°C
Power Supply:	13amp

MODEL	GROSS CAPACITY	EXTERNAL DIMENSIONS (H x W x D mm)	POWER CONSUMPTION	DECIBELS	WEIGHT (kg)	RRP
MED140	145L	819 x 598 x 648*	135W	53dBA	38	£749.00

* Depth 598mm without handle



OPTIONAL EXTRAS

M4-UCRSTACK	Stacking kit	£76.03
-------------	--------------	---------------

