


product information

Fully powered height adjustment 45-96cm	Fully powered multi-position head section +70° / -35°	Fully powered multi-position foot section +70° / -35°	Maximum user weight (35st.)	Easy to use handset
Retractable positioning castors	Adjustable height levelling foot	Fully upholstered head & foot sections	Double layer density foam for support & comfort	Luxury 65mm deep contoured upholstery
Anti-bacterial/anti-microbial vinyl	Paper roll holder (supplied as standard)	Easy clean design and materials	3 year warranty on all moving components	Lifetime warranty on frame

accessories * these accessories can only be fitted when ordered with the couch


Sun-FVHC/LEG1 Leg Supports	Sun-FVHC/PHB1 Phlebotomy Arm	Sun-FVHC/FSW1* Foot Switch	Sun-FVHC/FSW2* Foot Switch 2 Section	Sun-FVHC/FSW3* Foot Switch 3 Section
Sun-FVHC/SUP1 Gynae Stirrups	Sun-CHP1 Head Pillow	Sun-FVHC/CBH1* Breathing Hole	Sun-FVHC/CSR Folding Side Rails (pair)	Sun-FVHC/IV2 IV Pole Attachment
Sun-FVHC/RAIL Medi-Rail Attachment	Sun-INST1* Installation Service & Packaging Recycling	Sun-DISPOSAL1* Removal & Disposal of Old Couch		

dimensions



HEAVY ITEM

anti-bacterial vinyl upholstery colours - please specify when ordering please see inside back cover for full specification of our vinyl range



Fusion Classic Plus+ 2 & 3 Section Fully Powered Electric Couches

new

The new Fusion Classic Plus+ extends our already popular Fusion range, providing enhanced user functionality and easier patient positioning through its fully powered height, head and foot sections. Fusion Classic Plus+ is simply, strong, stable, aesthetically pleasing and easy to clean.

Fusion Classic Plus+

Fully powered height adjustment

Fully powered head section

Fully powered foot section

Code: **Sun-FVHC+3E**
Fusion Classic Plus+ 3 Section Electric Couch with powered height, head & foot section (Shown with optional accessories)

7-10 DAY DELIVERY*
Normally delivered within 7-10 working days (UK Only)



Code: **Sun-FVHC+2E**
Fusion Classic Plus+ 2 Section Electric Couch with powered height & head section