

**Combined Refrigerator
and Freezer**
HYCD-282A



**Refrigerator &
Freezer Combination**

Independent display for refrigerator
section and freezer section



Alarm Functions

High/low temperature, power failure,
senor error, low battery, high/low
ambient temperature, and remote
alarm



**Standard Stainless
Steel Interior**

Both refrigerator and freezer
chamber are stainless steel
interior, durable, easy to
clean and anti-corrosion

Ergonomic Design



Lockable
design



Intelligent
alarm



Foot/Caster
design



Test hole
design

Specifications

HYCD-282A

Technical Data

Type	Combined Refrigerator and Freezer
Cabinet Type	Upright
Climate Class	4
Cooling Type	Refrigerator: Forced air cooling Freezer: Direct cooling
Defrost Mode	Refrigerator: Auto Freezer: Manual
Refrigerant	CFC-Free
Noise Level (dB(A))	43

Performance

Temperature Range(°C)	Refrigerator: 2~8 Freezer: -20~-40
-----------------------	---------------------------------------

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220~240/50
Power(W)	400
Electrical Current(A)	2.8

Dimensions

Capacity(L/Cu.Ft)	Refrigerator:185/6.5	Freezer:97/3.4
Net/Gross Weight(approx)	145/160(kg)	319.7/352.7(lbs)
Interior Dimensions(W*D*H)	Refrigerator:605*510*720(mm) Freezer:515*465*440(mm)	23.8*20.0*28.3(in) 20.3*18.3*17.3(in)
Exterior Dimensions(W*D*H)	736*660*1810(mm)	29.0*26.0*71.3(in)
Packing Dimensions(W*D*H)	775*720*1930(mm)	30.5*28.3*76(in)
Container load (20'/40'/40'H)	21/45/45	



Functions

High/Low Temperature	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves/Drawers	Refrigerator: 3/1 Freezer: 0/2
USB Interface	Y

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R.China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.

