

## RLVF2025

### Basic Laboratory Freezer



With a robust steel exterior and easy clean plastic interior the large freezer from the Labcold Basic range is designed to provide laboratory storage at a lower price while still providing the confidence of the Labcold brand and support.

Lower cost production methods have been used in order to offer the best price possible without significant sacrifices in the important areas of function and reliability. This means your freezer still comes with a digital display, temperature alarms and a door lock as standard.



This roomy 543 litre freezer has six full width fixed shelves and swivel wheels to allow for easy manoeuvrability and ease of cleaning. A digital temperature display shows you the temperature inside the cabinet.

As you would expect from Labcold, it is backed by a comprehensive 1 year parts and labour warranty.\*

| Part Number                 | RLVF2025                                            |
|-----------------------------|-----------------------------------------------------|
| Capacity (Litres)           | 543                                                 |
| Style                       | Freestanding                                        |
| Exterior Dimensions (HxWxD) | 1960 x 750 x 755mm                                  |
| Interior Dimensions (HxWxD) | 1603 x 622 x 545mm                                  |
| Interior Construction       | ABS Plastic                                         |
| Exterior Construction       | Powder Coated Steel                                 |
| Net Weight                  | 100kg                                               |
| Energy Consumption          | 0.90kWh/24h                                         |
| Set Point                   | -20°C                                               |
| Temperature Range           | -18°C to -24°C                                      |
| Door Style                  | Solid                                               |
| Self Closing Door           | ✓                                                   |
| Door Hinging                | Right<br>(can be reversed at time of manufacture)   |
| Sparkfree                   | ✗                                                   |
| Access Port                 | Optional - order at time of manufacture             |
| Shelves                     | 6 fixed shelves<br>Usable dimensions<br>622 x 510mm |
| Wheels                      | 2 swivel and 2 swivel with brake                    |
| Door Lock                   | ✓                                                   |
| Manual Defrost              | ✓                                                   |
| Digital Temperature Display | ✓                                                   |
| High Temperature Alarm      | ✓                                                   |
| Low Temperature Alarm       | ✓                                                   |
| Door Open Alarm             | ✓                                                   |
| 1 Year Warranty*            | ✓                                                   |

\* labour warranty UK only.

01256 705 570

[www.labcold.com](http://www.labcold.com)

