

Medical & Mortuary Products



Product Name:

Hospital Pattern Disposal Sink - DU

Product Code:

VANTDU600x600RE

Product Description:

A press formed Stainless Steel Sluice with continuous flushing rim. No upstand, edges rimmed and turned down all round.
 Available with both Top Entry and Rear Entry feeds to the Sluice with or without Grid (See picture)
 Top Entry models supplied with cistern connecting tube.

Variations

VANTDU600X600TE - Hospital Pattern Disposal Sink – Top Entry
 BSGRID - Disposal Sink Grid and Mounting Pillar Assembly

Technical Specification:

<i>Construction:</i>	Pressed, Folded and welded.	<i>Net Weight</i>	6.5kg
<i>Material:</i>	304 Stainless Steel	<i>Sluice Outlet:</i>	110mm DIA
<i>Material Thickness:</i>	1.2mm	<i>Fixings:</i>	Canti-lever Wall Brackets

