

Medical & Mortuary Products



Product Name:

Hospital Pattern Disposal Unit - DUHS

Product Code:

VANTDUHS1600x600RELH

Product Description:

A press formed plain top of Stainless Steel incorporating a rectangular sink bowl without overflow or tap holes and a Sluice with continuous flushing rim. No upstand, edges rimmed and turned down all round.

Available with both Top Entry and Rear Entry feeds to the Sluice and Right Hand or Left Hand Sluice Unit position.

Top Entry models supplied with cistern connecting tube.

Variations

VANTDUHS1600x600TERH – Top Entry, Right Hand Sluice
 VANTDUHS1600x600RERH – Rear Entry, Right Hand Sluice
 VANTDUHS1600x600TELH – Top Entry, Left Hand Sluice

Technical Specification:

Net Weight:	15.5kg	Bowl Size:	500mm x 400mm
Material:	304 Stainless Steel	Sluice Outlet:	110mm DIA
Material Thickness:	1.2mm	Fixings:	Canti-lever Wall Bracket

